

MedRound Infinity Cart

Minimise errors. Enhance workflow. Reduce bending.



Dispense accurately, efficiently and comfortably with the new high-capacity, height adjustable MedRound.

Store More

- Enhance your workflow and reduce the need for restocking with up to 8 drawers of visible, organised and secure storage.
- A centralised bottle storage drawer pulls out to allow bottles to be contained separately.
- Drawers configured to suit you - with different depth options available.
- Medication is securely housed and protected by high-security Camlock, Smartlock or Digilock systems.

Bend Less

Long medication round solutions have traditionally been fixed height products with drawers starting from the base up.

We have worked with nurses to design a unique height adjustable solution which brings every dispensing tray within easy reach, without the need for bending.

24/7 Access to ePMA

- All-in-one touchscreen computer with high-speed processor allows instant access to ePMA at the point of care.
- Dual hot-swappable batteries allow the workstation to be in constant use. Batteries can be replaced in seconds.
- The battery management system provides real time information to the user and IT colleagues.

Perfect for Infection Control

- IP65 front panel protection and fanless All-in-One computer together with smooth contours on the cart make this the perfect solution for infection control.

Accessories



Scanner



Mobile Printer



Sharps Bin



Infection Control
Keyboard and Mouse



Drawer Inserts



Made in the UK



Project management



Phased delivery



Live chat service





Preimaged hard drives



Device configuration

MedRound Infinity Cart

Minimise errors. Enhance workflow. Reduce bending.




Cart

Work Surface	958 X 478 mm
Manoeuvring	Contoured handles and 4 large, lockable dual wheel castors ensure the cart is easy to manoeuvre.
Height adjustment stroke	354mm
Storage capacity	Up to 78 Litres, depending on drawer config
Storage Security	Secure metal enclosure, protected by your choice of lock.
Drawer Dimension Options	427 L x 312 W x 75, 150 or 225 Depth (mm)

Lock Configurations


HIGH-SECURITY CAMLOCK



DIGILOCK



SMARTLOCK



Medically Certified

EN60601	Certified Power Supply
Infection Control	No external fixings No visible vents

All-in-One Computer

Size	19", 22" or 24"
Resolution	1280x1080 (19") or 1920x1080 (22"/24")
View Angle	170°(H)/160°(V) or 178°(H)/178°(V)
Processor	Intel® Quad-Core J1900 / Intel® Skylake Dual-Core i5
OS Support	Windows 7, 8.1 and 10
Power Requirement	AC100 V~240V. DC24V
Operating Temperature	0°C ~ 40°C (32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C (-4°F ~ 140°F)
Storage Humidity	10%~95% @ 35°, non-condensing
Degree of Protection	Front Panel IP-65 Certified, IP x 1
System Cooling	Fanless
Net Weight	8KG

Approved Cleaning Solutions

Made with materials that are resistant to cleaning solutions, the Cart can be cleaned with most hospital-standard cleaning solutions. A list of approved cleaning solutions is available on request. Please contact us for a full list