

# VersaLite ePMA Infinity Cart

## Revolutionise Medicine Dispensing



The new VersaLite ePMA Infinity Cart has been developed to provide the ultimate combination of functionality, mobility and user comfort.

Improve workflow and patient safety with timely and accurate electronic prescribing.

### 24/7 access to patient records

- All-in-one touchscreen computer with Intel® SkyLake dual core i5 processor allows instant access to patient records at the point of care.
- Hot "swappable" batteries allow the cart to be in constant use.
- The battery management system provides real-time information to the user and IT colleagues. Batteries can be replaced in seconds.

### Lightweight and compact design

- Small footprint makes the cart easy to position in busy wards.
- Conveniently placed grab handles allow precise and easy movement.
- Large rear castors provide smooth and stable movement even across uneven surfaces such as thresholds.
- Designed for infection control and ease of cleaning.
- Digilock or secure camlock ensures secure storage of medication.
- Drawer configuration can be adapted to suit specific workflow. Up to 32 litres of storage capacity.

### Designed for individuals

- Simple hand operated height adjustment.
- Ergonomically designed base allows the user to find their perfect standing or seated position.
- Retractable keyboard tray features a simple lift-and-release mechanism allowing the user to set the height, distance and angle independently of the workstation, so they maintain a healthy posture.

## Accessories



Scanner




Sharps bin



Mobile Printer

# VersaLite ePMA Infinity Cart

## Revolutionise Medicine Dispensing



### Cart

<b>Work Surface</b>	513 X 636mm
<b>Manoeuvring</b>	Contoured handles and larger rear for ease of manoeuvre.
<b>Keyboard</b>	Height and angle adjustable. 180 degree lateral movement
<b>Wheels</b>	2x front wheels with brakes, 2x larger rear wheels for additional stability
<b>Height adjustment stroke</b>	354mm

### Medically Certified


<b>EN60601</b>	Certified Power Supply
<b>Infection Control</b>	No External fittings No Visible vents

### Display

<b>Size</b>	19", 22" or 24"
<b>Resolution</b>	1900x1080
<b>View Angle</b>	170°(H)/160°(V) or 178°(H)/178°(V)
<b>Processor</b>	Intel® Skylake Dual-Core i5
<b>System Memory</b>	DDR4 up to 32GB

### Mechanical and Environmental

<b>Power Requirement</b>	AC100 V~240V. DC24V
<b>Operating Temperature</b>	0°C ~ 40°C (32°F ~ 104°F)
<b>Storage Temperature</b>	-20°C ~ 60°C (-4°F ~ 140°F)
<b>Storage Humidity</b>	10%~95% @ 35°, non-condensing
<b>Degree of Protection</b>	Front Panel IP-65 Certified, IP x 1
<b>System Cooling</b>	Fan-less
<b>Net Weight</b>	8KG
<b>EMC</b>	CE/FCC Class B, UL-60601-1, EN-60601-1



### Approved Cleaning Solutions

- Quaternary ammonium chloride compounds (such as AirX109)
- Bleach-Household strength (1:10 minimum dilution strength)
- Ethyl or Isopropyl alcohol-70%
- Hydrogen Peroxide
- Phenolic disinfectants (such as Lyol, WexCide)
- Ovation
- Fantastic
- Formula 409
- 1.6% Aqueous Ammonia "Green Soap" USP (United States Pharmacopoeia)
- Clinell Wipes
- Chlor Clea