

medDispense® CS series

Expanded capacity for bulk supply



Overview

The medDispense® CS series automated dispensing cabinets provide you with expanded capacity to secure bulk supply or other large items.

- ✓ Capacity to hold from 600 to 1,200 items
- ✓ Combine with medGuard™ software to securely store, trace and dispense controlled substances
- ✓ Powered by medLogic™ to deliver increased efficiency and flexibility for your team and a safer environment for your patients

Features

- Provides expanded capacity to secure bulk supply or other large items
- Combine with the medGuard™ software module for secure controlled substance storage and closed-loop medication restocking
- Integrates with C and F series medDispense® automated dispensing cabinets
- Choose from clear front or steel supply doors
- Microban® antimicrobial infused surfaces
- Comprehensive 5- year warranty

Other unit sizes in this series



Recommended for

- Psychiatric
- Critical access hospitals
- Simulation labs
- Long-term acute care

Compatible with



C series



F series

medDispense® CS series

Expanded capacity for bulk supply

Accessories



Bulk storage bins



Overbridge

Specifications

Max unit capacity	600 or 1,200
Width	686 mm
Depth	559 mm (compact unit) or 711 mm (full depth unit)
Height	1,118 mm (600-capacity unit) or 1,981 mm (1,200-capacity unit)
Weight	74-85 kgs (600-capacity) or 120-136 kgs (1,200-capacity unit)
Max number of drawers	12, 30, or 60 (dependent on unit capacity)
Material construction	Steel frame with polymer or aluminium doors (1,200-unit capacity option only)
Power	120VAC @ 60Hz

Software

- All medDispense® automated dispensing cabinets utilize the medLogic™ software platform
- The medLogic™ software platform delivers increased efficiency and flexibility for your team and a safer environment for your patients.
- **Efficiency:**
 - Single medication dispensing across multiple medDispense® automated dispensing cabinets.
 - Restock faster with barcode stocking.
 - FIFO stocking workflow ensures oldest stocked items are dispensed first, which helps reduced expired medication expenses.
- **Safety:**
 - Fully supports the UK NHS ECD register.
 - On-demand stock report generation by item/medication.
 - Help reduce medication errors with our configurable medication search.
 - Increase medication traceability with lot number stocking.
- **Flexibility:**
 - Fractional dose stocking & dispensing helps reduce medication waste.
 - High definition 16:9 and laptop display with the PharmacyCenter™ and NurSync™ modules provides increased technology flexibility for the pharmacy team.
- medLogic™ secures and tracks both your medications and supplies while seamlessly integrating with any hospital information system.
- medLogic™ modules improve workflow efficiency, safety and visibility so you can save time and focus on delivering exceptional patient care.
- **NurSync™ Remote Medication Queuing:**
 - Enables nurses to create medication queues for their patients to help reduce medication retrieval times, create a more efficient workflow and spend more time with patients.
- **medGuard™ Controlled Drug Inventory Management:**
 - Helps close the loop of controlled drug medication delivery and monitors for potential diversion with controlled restocking access and software safeguards.
 - Discrepancy tracking ensures the retrieval quantity recorded in medGuard™ matches the quantity restocked in the medDispense® automated dispensing cabinet.
- **PharmacyCenter™ Reporting and System Administration:**
 - Enables pharmacy staff to manage their medications across multiple medDispense® automated dispensing cabinets from one easy-to-use dashboard and reporting system.
 - Manage and maintain your facility's automated dispensing cabinets from any PC.
 - Visibility and management of your inventory across all automated dispensing cabinets.
 - Report on key metrics that matter to your facility, such as: inventory, top usage, transactions, controlled medications, and more.

