

# medDispense® F series

The highest level of security with single dose dispensing

## Overview

The medDispense® F series automated dispensing cabinets' single dose drawer technology allows medical staff to access only the specified compartment within the drawer for maximum security of items like controlled drugs.

- ✓ Locate medication quickly with illuminated drawers that light up when medication is ready to dispense
- ✓ Available in multiple configurations with up to 15 drawer sizes to accommodate all types of medications
- ✓ The medLogic™ software platform delivers increased efficiency and flexibility for your pharmacy, caregivers, and IT team
- ✓ Countertop model provides all the features of the larger unit in a smaller footprint to accommodate your medication security needs with less space
- ✓ With a maximum capacity of 384 items in 8 customizable levels, the countertop model offers plenty of options for facilities that might be short on space—and budget
- ✓ Includes an electronic CD register, designed for the NHS

## Features

- Provides the highest level of security with single-dose steel progressive drawers
- Due to durability, F series medDispense® automated dispensing cabinets are approved for use in HM Prison Service
- Pick-to-Light guidance technology guides you to the drawer with the selected medication
- Steel construction fully supports the NHS controlled drug requirements
- Supports three stocking modes: single, mixed, and linear
- Medication management in all the sizes you need
- Manufactured in the UK
- Comprehensive 5-year warranty

## Recommended for

- Hospitals
- Ambulatory surgery centres
- Long-term acute care



## Other unit sizes in this series



## Compatible with



C series



CS series



CT series



V series

# medDispense® F series

The highest level of security with single dose dispensing

## Accessories



Remote lock



Remote lock with temperature monitoring



External return bin



Barcode scanner



Thermal printer



Label printer



Fingerprint scanner

## Specifications

Max unit capacity	1,248
Width	545 mm
Depth	520 mm
Height	2,003 mm
Weight	245 kgs
Material construction	Steel frame and drawers
Power	120VAC @ 60Hz

## Software

- All medDispense® automated dispensing cabinets utilize the medLogic™ software platform
- The medLogic™ software platform delivers increased efficiency and flexibility for your team and a safer environment for your patients.
- **Efficiency:**
  - Single medication dispensing across multiple medDispense® automated dispensing cabinets.
  - Restock faster with barcode stocking.
  - FIFO stocking workflow ensures oldest stocked items are dispensed first, which helps reduced expired medication expenses.
- **Safety:**
  - Fully supports the UK NHS ECD register.
  - On-demand stock report generation by item/medication.
  - Help reduce medication errors with our configurable medication search.
  - Increase medication traceability with lot number stocking.
- **Flexibility:**
  - Fractional dose stocking & dispensing helps reduce medication waste.
  - High definition 16:9 and laptop display with the PharmacyCenter™ and NurSync™ modules provides increased technology flexibility for the pharmacy team.
- medLogic™ secures and tracks both your medications and supplies while seamlessly integrating with any hospital information system.
- medLogic™ modules improve workflow efficiency, safety and visibility so you can save time and focus on delivering exceptional patient care.
- **NurSync™ Remote Medication Queuing:**
  - Enables nurses to create medication queues for their patients to help reduce medication retrieval times, create a more efficient workflow and spend more time with patients.
- **medGuard™ Controlled Drug Inventory Management:**
  - Helps close the loop of controlled drug medication delivery and monitors for potential diversion with controlled restocking access and software safeguards.
  - Discrepancy tracking ensures the retrieval quantity recorded in medGuard™ matches the quantity restocked in the medDispense® automated dispensing cabinet.
- **PharmacyCenter™ Reporting and System Administration:**
  - Enables pharmacy staff to manage their medications across multiple medDispense® automated dispensing cabinets from one easy-to-use dashboard and reporting system.
  - Manage and maintain your facility's automated dispensing cabinets from any PC.
  - Visibility and management of your inventory across all automated dispensing cabinets.
  - Report on key metrics that matter to your facility, such as: inventory, top usage, transactions, controlled medications, and more.

