

medDispense® V series

Storage and traceability of bulk medications and supplies



Overview

The medDispense® V series automated dispensing cabinets provide storage and traceability of bulk medications and supplies. Ideal for central pharmacy controlled drug storage.

- ✓ Easily configure cabinets for your needs with three standard door sizes that can be mixed within one frame
- ✓ Locate medication quickly with illuminated cabinet and doors
- ✓ Enhance visibility and reduce waste with clear door options
- ✓ Capacity to hold up to 720 items
- ✓ Combine with medGuard™ software to securely store, trace and dispense controlled substances
- ✓ Powered by medLogic™ to deliver increased efficiency and flexibility for your team and a safer environment for your patients

Features

- Provides storage and traceability of bulk medications and supplies
- Ideal for central pharmacy controlled drug storage
- Combined with the medGuard™ module, the V series provides a closed loop solution for your controlled drug management
- Integrates with C and F series medDispense® automated dispensing cabinets
- Items are easily located by pop open doors or guiding lights for quick and intuitive access
- Microban® antimicrobial infused surfaces
- Comprehensive 5-year warranty

Recommended for

- Hospitals
- Rehabilitation
- Long-term acute care
- Ambulatory surgery centres

Other unit sizes in this series



Compatible with



C series



F series

medDispense® V series

Storage and traceability of bulk medications and supplies

Specifications

| | |
|-----------------------|--|
| Max unit capacity | 720 |
| Width | 1,093 mm |
| Depth | 500 mm |
| Height | 2,030 mm |
| Weight | 55 kgs (single locker) or 175 kgs (triple locker) |
| Material construction | Steel frame with either clear polycarbonate or steel doors |
| Power | 120VAC @ 60Hz |

Software

- All medDispense® automated dispensing cabinets utilize the medLogic™ software platform
- The medLogic™ software platform delivers increased efficiency and flexibility for your team and a safer environment for your patients.
- **Efficiency:**
 - Single medication dispensing across multiple medDispense® automated dispensing cabinets.
 - Restock faster with barcode stocking.
 - FIFO stocking workflow ensures oldest stocked items are dispensed first, which helps reduced expired medication expenses.
- **Safety:**
 - Fully supports the UK NHS ECD register.
 - On-demand stock report generation by item/medication.
 - Help reduce medication errors with our configurable medication search.
 - Increase medication traceability with lot number stocking.
- **Flexibility:**
 - Fractional dose stocking & dispensing helps reduce medication waste.
 - High definition 16:9 and laptop display with the PharmacyCenter™ and NurSync™ modules provides increased technology flexibility for the pharmacy team.
- medLogic™ secures and tracks both your medications and supplies while seamlessly integrating with any hospital information system.
- medLogic™ modules improve workflow efficiency, safety and visibility so you can save time and focus on delivering exceptional patient care.
- **NurSync™ Remote Medication Queuing:**
 - Enables nurses to create medication queues for their patients to help reduce medication retrieval times, create a more efficient workflow and spend more time with patients.
- **medGuard™ Controlled Drug Inventory Management:**
 - Helps close the loop of controlled drug medication delivery and monitors for potential diversion with controlled restocking access and software safeguards.
 - Discrepancy tracking ensures the retrieval quantity recorded in medGuard™ matches the quantity restocked in the medDispense® automated dispensing cabinet.
- **PharmacyCenter™ Reporting and System Administration:**
 - Enables pharmacy staff to manage their medications across multiple medDispense® automated dispensing cabinets from one easy-to-use dashboard and reporting system.
 - Manage and maintain your facility's automated dispensing cabinets from any PC.
 - Visibility and management of your inventory across all automated dispensing cabinets.
 - Report on key metrics that matter to your facility, such as: inventory, top usage, transactions, controlled medications, and more.

